

Product datasheet

Anti-KIF26A antibody ab57299

1 图像

概述

| | |
|-------|---|
| 产品名称 | Anti-KIF26A抗体 |
| 描述 | 小鼠单克隆抗体to KIF26A |
| 经测试应用 | 适用于: WB |
| 种属反应性 | 与反应: Human |
| 免疫原 | Recombinant full length protein, corresponding to amino acids 1-204 of Human KIF26A |

性能

| | |
|------|---|
| 形式 | Liquid |
| 存放说明 | Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. |
| 存储溶液 | Preservative: None PBS, pH 7.2 |
| 纯度 | Protein G purified |
| 克隆 | 单克隆 |
| 同种型 | IgG1 |
| 轻链类型 | kappa |

应用

Our [Abpromise guarantee](#) covers the use of **ab57299** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| 应用 | Ab评论 | 说明 |
|----|------|----|
|----|------|----|

WB


应用说明 WB: Use at a concentration of 1-5 µg/ml.

Not yet tested in other applications.
Optimal dilutions/concentrations should be determined by the end user.

靶标

| | |
|--------------|---|
| 功能 | Atypical kinesin that plays a key role in enteric neuron development. Acts by repressing a cell growth signaling pathway in the enteric nervous system development, possibly via its interaction with GRB2 that prevents GRB2-binding to SHC, thereby attenuating the GDNF-Ret signaling. Binds to microtubules but lacks microtubule-based motility due to the absence of ATPase activity. |
| 序列相似性 | Belongs to the kinesin-like protein family. KIF26 subfamily. Contains 1 kinesin-motor domain. |
| 细胞定位 | Cytoplasm > cytoskeleton. |

图片



Predicted band size : 23 kDa
KIF26A antibody (ab57299) at 1ug/lane +
HeLa cell lysate at 25ug/lane.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors