

## Product datasheet

# Rabbit F(ab')<sub>2</sub> Anti-Goat IgG H&L (FITC) preadsorbed ab5750

### 概述

<b>产品名称</b>	兔F(ab') <sub>2</sub> 抗山羊IgG H&L (FITC)预吸附二抗
<b>宿主</b>	Rabbit
<b>靶标种属</b>	Goat
<b>特异性</b>	Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Fluorescein, anti-Rabbit Serum, Goat IgG and Goat Serum. No reaction was observed against anti-Pepsin, anti-Rabbit IgG F(c) or Human Serum Proteins.
<b>经测试应用</b>	<b>适用于:</b> ICC/IF, Immunomicroscopy, Flow Cyt, ELISA
<b>最小交叉反应</b>	Human <a href="#">more details</a>
<b>免疫原</b>	Goat IgG whole molecule.
<b>偶联物</b>	FITC. Ex: 493nm, Em: 528nm

### 性能

<b>形式</b>	Liquid
<b>存放说明</b>	Shipped at 4°C. Store at +4°C.
<b>存储溶液</b>	Preservative: 0.01% Sodium Azide Constituents: 10mg/ml BSA (IgG and protease free), 0.15M Sodium Chloride, 0.02M Potassium Phosphate. pH 7.2
<b>纯化说明</b>	This product was prepared from monospecific antiserum by immunoaffinity chromatography using Goat IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities, pepsin digestion and chromatographic separation.
<b>共轭说明</b>	Fluorescein isothiocyanate (FITC) (Molecular Weight 390 daltons) Absorption Wavelength: 495 nm Emission Wavelength: 528 nm Fluorochrome/Protein Ratio: 3.3 moles FITC per mole of Rabbit IgG F(ab') <sub>2</sub>
<b>克隆</b>	多克隆
<b>同种型</b>	IgG

### 应用

Our [Abpromise guarantee](#) covers the use of **ab5750** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
ICC/IF		1/1000 - 1/5000.
Immunomicroscopy		Use at an assay dependent dilution.
Flow Cyt		1/500 - 1/2500.
ELISA		1/10000 - 1/50000.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors