

Product datasheet

Anti-Protocadherin beta 3 antibody ab57598

1 图像

概述

产品名称	Anti-Protocadherin beta 3抗体
描述	小鼠单克隆抗体to Protocadherin beta 3
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	Recombinant fragment: ASEEIRKTFQ LNPITGDMQL VKYLNFEAIN SYEVDIEAKD GGGLSGKSTV IVQVVDVNDN PPELTLSSVN SPIPENSGET VLAVFSVSDL DSGDNGRV, corresponding to amino acids 284-382 of Human Protocadherin beta 3 <a href="#">Run BLAST with ExPASy</a> <a href="#">Run BLAST with NCBI</a>

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None PBS, pH 7.2
纯度	Protein G purified
克隆	单克隆
同种型	IgG2a
轻链类型	kappa

应用

Our [Abpromise guarantee](#) covers the use of **ab57598** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
----	------	----

WB

应用说明 WB: Use at a concentration of 1-5 µg/ml.

This antibody has only been tested in WB against the recombinant fragment used as immunogen. We have no data on the detection of endogenous protein.


Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

## 靶标

功能	Potential calcium-dependent cell-adhesion protein. May be involved in the establishment and maintenance of specific neuronal connections in the brain.
序列相似性	Contains 6 cadherin domains.
细胞定位	Cell membrane.

## 图片



Western blot against tagged recombinant protein immunogen using ab57598  
Protocadherin beta 3 antibody at 1ug/ml.  
Predicted band size of immunogen is 37 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors