

Product datasheet

Anti-PDCL3 antibody ab57708

1 图像

概述

产品名称	Anti-PDCL3抗体
描述	小鼠单克隆抗体to PDCL3
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	Recombinant full length protein, corresponding to amino acids 1-240 of Human PDCL3

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None PBS, pH 7.2
纯度	Protein G purified
克隆	单克隆
同种型	IgG1
轻链类型	kappa

应用

Our [Abpromise guarantee](#) covers the use of **ab57708** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
----	------	----

WB


应用说明 WB: Use at a concentration of 1-5 µg/ml.

Not yet tested in other applications.
Optimal dilutions/concentrations should be determined by the end user.

靶标

功能	Modulates the activation of caspases during apoptosis. Is a substrate for <i>Orgyia pseudotsugata</i> multicapsid polyhedrosis virus (OpMNPV) IAP-mediated ubiquitination.
组织特异性	Expressed in all tissues examined including spleen, thymus, prostate, testis, ovary, small intestine and colon.
序列相似性	Belongs to the phosphatidylethanolamine phosphatase family.
细胞定位	Cytoplasm.

图片



Predicted band size : 29 kDa

PDCL3 antibody (ab57708) at 1ug/lane +
Jurkat cell lysate at 25ug/lane.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors