

Product datasheet

Anti-TSC22 domain family, member 4 antibody ab57885

1 图像

概述

|       |  |
|-------|--|
| 产品名称  | Anti-TSC22 domain family, member 4抗体   |
| 描述    | 小鼠单克隆抗体到 TSC22 domain family, member 4   |
| 宿主    | Mouse  |
| 经测试应用 | 适用于: WB  |
| 种属反应性 | 与反应: Human   |
| 免疫原   | Recombinant full length protein, corresponding to amino acids 1-396 of Human TSC22 domain family, member 4 |

性能

|      |   |
|------|---|
| 形式   | Liquid  |
| 存放说明 | Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. |
| 存储溶液 | Preservative: None<br>PBS, pH 7.2   |
| 纯度   | Protein G purified  |
| 克隆   | 单克隆   |
| 同种型  | IgG1  |
| 轻链类型 | kappa   |

应用

Our [Abpromise guarantee](#) covers the use of **ab57885** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| 应用 | Ab评论 | 说明  |
|----|------|---|
| WB |      | Use a concentration of 1 - 5 µg/ml. Predicted molecular weight: 41 kDa. |


靶标

|       |                                       |
|-------|---------------------------------------|
| 功能    | Transcriptional repressor.            |
| 序列相似性 | Belongs to the TSC-22/Dip/Bun family. |
| 细胞定位  | Nucleus.                              |

---

## 图片

---



TSC22 domain family member 4 antibody (ab57885) at 1ug/lane + K-562 cell lysate at 25ug/lane.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors