

Product datasheet

Anti-DNER antibody ab57949

1 图像

概述

产品名称	Anti-DNER抗体
描述	小鼠单克隆抗体to DNER
经测试应用	适用于: WB
种属反应性	与反应: Rat, Human
免疫原	Recombinant fragment, corresponding to amino acids 368-477 of Human DNER
阳性对照	PC-12 cell lysate

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None PBS, pH 7.2
纯度	Protein G purified
克隆	单克隆
同种型	IgG3
轻链类型	kappa

应用

Our [Abpromise guarantee](#) covers the use of **ab57949** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		Use a concentration of 1 - 5 µg/ml. Detects a band of approximately 90-100 kDa (predicted molecular weight: 78.4 kDa). Larger band size is due to glycosylation.

靶标

功能	Activator of the NOTCH1 pathway. May mediate neuron-glia interaction during astrocytogenesis.
组织特异性	Expressed in brain, spinal cord and adrenal gland.
序列相似性	Contains 10 EGF-like domains. Contains 1 follistatin-like domain.
细胞定位	Cell membrane. Present on the membrane of dendrites and cell bodies but excluded from axonal membrane. Also found in early endosomes in the somatodendritic region (By similarity).

Anti-DNER antibody 图像



Predicted band size : 78.4 kDa

Western Blot analysis of DNER antibody (ab57949).

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors