

Product datasheet

Anti-IBRDC1 antibody ab58041

1 图像

概述

产品名称	Anti-IBRDC1抗体
描述	小鼠单克隆抗体to IBRDC1
宿主	Mouse
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	Recombinant fragment: KVQLGQVEIK CPITECFEFL EETTVVYNLT HEDSIKYKYF LELGRIDSST KPCPQCKHFT TFKKKGHIPT PSRSESKYKI QCPTCQFVWC FKCHSPWHE, corresponding to amino acids 2-101 of Human IBRDC1 <a href="#">Run BLAST with ExPASy</a> <a href="#">Run BLAST with NCBI</a>

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None PBS, pH 7.2
纯度	Protein G purified
克隆	单克隆
同种型	IgG2a
轻链类型	kappa

应用

Our [Abpromise guarantee](#) covers the use of **ab58041** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		Use a concentration of 1 - 5 µg/ml. Predicted molecular weight: 32 kDa.

## 靶标

相关性	IBRDC1 is an E3 ubiquitin protein ligase which accepts ubiquitin from E2 ubiquitin conjugating enzymes in the form of a thioester and then directly transfers the ubiquitin to targeted substrates.
细胞定位	Membrane; Single pass membrane protein.

## 图片



IBRDC1 antibody (ab58041) at 1ug/lane +  
HepG2 cell lysate at 25ug/lane.

Western blot - IBRDC1 antibody (ab58041)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors