

Product datasheet

Anti-IBRDC1 antibody ab58041

1 图像

概述

| | |
|-------|---|
| 产品名称 | Anti-IBRDC1抗体 |
| 描述 | 小鼠单克隆抗体to IBRDC1 |
| 宿主 | Mouse |
| 经测试应用 | 适用于: WB |
| 种属反应性 | 与反应: Human |
| 免疫原 | Recombinant fragment: KVQLGQVEIK CPITECFEFL EETTVVYNLT HEDSIKYKYF LELGRIDSST KPCPQCKHFT TFKKKGHIPT PSRSESKYKI QCPTCQFVWC FKCHSPWHE, corresponding to amino acids 2-101 of Human IBRDC1 Run BLAST with ExPASy Run BLAST with NCBI |

性能

| | |
|------|---|
| 形式 | Liquid |
| 存放说明 | Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. |
| 存储溶液 | Preservative: None PBS, pH 7.2 |
| 纯度 | Protein G purified |
| 克隆 | 单克隆 |
| 同种型 | IgG2a |
| 轻链类型 | kappa |

应用


Our [Abpromise guarantee](#) covers the use of **ab58041** in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| 应用 | Ab评论 | 说明 |
|----|------|---|
| WB | | Use a concentration of 1 - 5 µg/ml. Predicted molecular weight: 32 kDa. |

靶标

| | |
|------|---|
| 相关性 | IBRDC1 is an E3 ubiquitin protein ligase which accepts ubiquitin from E2 ubiquitin conjugating enzymes in the form of a thioester and then directly transfers the ubiquitin to targeted substrates. |
| 细胞定位 | Membrane; Single pass membrane protein. |

图片



IBRDC1 antibody (ab58041) at 1ug/lane +
HepG2 cell lysate at 25ug/lane.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors