

Product datasheet

Anti-OAS2 antibody ab60091

1 图像

概述

产品名称	Anti-OAS2抗体
描述	兔多克隆抗体to OAS2
宿主	Rabbit
经测试应用	适用于: ELISA, IHC-P
种属反应性	与反应: Human
免疫原	Synthetic peptide (10-30 aa in length) at the C-term of last 50 aa of human OAS2 conjugated to KLH
阳性对照	Human breast carcinoma tissue.

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
存储溶液	Preservative: 0.09% Sodium azide Constituent: PBS
纯度	Protein G purified
纯化说明	ab60091 is purified through a protein G column and eluted out with both high and low pH buffers and neutralized immediately after elution then followed by dialysis against PBS.
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee

Abpromise™承诺保证使用ab60091于以下的经测试应用

“应用说明”部分下显示的仅为推荐的起始稀释度；实际最佳的稀释度/浓度应由使用者检定。


应用	Ab评论	说明
ELISA		1/1000.

应用	Ab评论	说明
IHC-P		1/50 - 1/100.

靶标

功能	May play a role in mediating resistance to virus infection, control of cell growth, differentiation, and apoptosis.
序列相似性	Belongs to the 2-5A synthase family.
细胞定位	Mitochondrion. Nucleus. Microsome. Endoplasmic reticulum. Associated with different subcellular fractions such as mitochondrial, nuclear, and rough/smooth microsomal fractions.

图片



Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors