

## Product datasheet

# Anti-CD90 / Thy1 antibody [30-H12] (FITC) ab62009

### 概述

产品名称	Anti-CD90 / Thy1 抗体[30-H12] (FITC)
描述	大鼠单克隆抗体[30-H12] to CD90 / Thy1 (FITC)
宿主	Rat
偶联物	FITC. Ex: 493nm, Em: 528nm
经测试应用	适用于: WB, IP, Flow Cyt, IHC-Fr
种属反应性	与反应: Mouse
免疫原	Full length mouse CD90/Thy1 protein.

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Preservative: 0.1% Sodium Azide Constituents: 40% Glycerol, 1% BSA, PBS, pH 7.2
纯度	Protein G purified
纯化说明	ab62009 was purified from tissue culture supernatant.
克隆	单克隆
克隆编号	30-H12
同种型	IgG2b

### 应用

Our [Abpromise guarantee](#) covers the use of **ab62009** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab 评论	说明
WB		Use at an assay dependent concentration. Predicted molecular weight: 18 kDa. This product is also suitable for use in identification of mouse T-cells subset population and phenotype analysis as a marker for progenitor cells during developmental stage.

应用	Ab 评论	说明
----	-------	----

IP Use at an assay dependent concentration.

Flow Cyt Use at an assay dependent concentration.  
(cell sorting or panning)

[ab136125](#) - Rat monoclonal IgG2b, is suitable for use as an isotype control with this antibody.

IHC-Fr Use at an assay dependent concentration.

## 靶标

功能	May play a role in cell-cell or cell-ligand interactions during synaptogenesis and other events in the brain.
序列相似性	Contains 1 Ig-like V-type (immunoglobulin-like) domain.
细胞定位	Cell membrane.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors