


Product datasheet

Anti-HOXB9 antibody ab65063

1 图像

概述

产品名称	Anti-HOXB9抗体
描述	兔多克隆抗体to HOXB9
经测试应用	适用于: WB, ELISA
种属反应性	与反应: Human 预测可用于: Mouse 
免疫原	Synthetic peptide (Human) from an internal sequence.
阳性对照	Extracts from K562 cells

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at -20°C. Stable for 12 months at -20°C.
存储溶液	Preservative: 0.02% Sodium Azide Constituents: 50% Glycerol, PBS (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), 150mM Sodium chloride, pH 7.4
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab65063** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/500 - 1/1000. Detects a band of approximately 28 kDa (predicted molecular weight: 28 kDa).
ELISA		1/10000.

靶标

<b>功能</b>	Sequence-specific transcription factor which is part of a developmental regulatory system that provides cells with specific positional identities on the anterior-posterior axis.
<b>序列相似性</b>	Belongs to the Abd-B homeobox family. Contains 1 homeobox DNA-binding domain.
<b>发展阶段</b>	Expressed in whole embryos and fetuses at 5-9 weeks from conception.
<b>细胞定位</b>	Nucleus.

## 图片



**All lanes** : Anti-HOXB9 antibody (ab65063)  
at 1/500 dilution

**Lane 1** : Extracts from K562 cells with no immunising peptide

**Lane 2** : Extracts from K562 cells with immunising peptide at 5 µg

Lysates/proteins at 5 µg per lane.

**Predicted band size** : 28 kDa

**Observed band size** : 28 kDa

**Additional bands at** : >117 kDa. We are unsure as to the identity of these extra bands.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors