

Product datasheet

Anti-A1BG antibody ab66011

1 图像

概述

产品名称	Anti-A1BG抗体
描述	兔多克隆抗体to A1BG
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	Synthetic peptide: VTFLLRREGD HEFLEVPEAQ EDVEATFPVH QPGNYSCSYR TDGEGALSEP , corresponding to amino acids 146-195 of Human A1BG Run BLAST with ExPASy Run BLAST with NCBI
阳性对照	HepG2 cell lysate.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Preservative: None Constituents: 2% Sucrose, PBS
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab66011** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


应用	Ab 评论	说明
----	-------	----

WB	Use a concentration of 1.25 µg/ml. Detects a band of approximately 54 kDa (predicted molecular weight: 54 kDa). Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.
----	--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

靶标

组织特异性	Plasma.
序列相似性	Contains 5 Ig-like V-type (immunoglobulin-like) domains.
细胞定位	Secreted.

图片



Anti-A1BG antibody (ab66011) at 1.25 µg/ml
+ HepG2 cell lysate at 10 µg

Secondary

HRP conjugated anti-Rabbit IgG at 1/50000
dilution

Predicted band size : 54 kDa

Observed band size : 54 kDa

Additional bands at : 30 kDa (possible non-specific binding).

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors