

Product datasheet

Anti-SLAP2 antibody ab66060

1 图像

概述

产品名称	Anti-SLAP2抗体
描述	小鼠单克隆抗体to SLAP2
宿主	Mouse
经测试应用	适用于: ELISA, WB
种属反应性	与反应: Human
免疫原	Recombinant full length human SLAP2 (AAH42041, 1 a.a. ~ 262 a.a) with tag.
阳性对照	SLAP2 transfected 293T cell lysate.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: None Constituents: 1X PBS, pH 7.2
纯度	Protein A purified
克隆	单克隆
同种型	IgG1
轻链类型	kappa

应用

Our [Abpromise guarantee](#) covers the use of **ab66060** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


应用	Ab评论	说明
ELISA		Use at an assay dependent dilution. Detection limit for recombinant tagged SLAP2 is approximately 0.03ng/ml if used as a capture antibody.
WB		Use a concentration of 1 - 5 µg/ml. Detects a band of approximately 29 kDa (predicted molecular weight: 29 kDa).

## 靶标

---

## 图片

---



**All lanes** : Anti-SLAP2 antibody (ab66060) at 1 µg/ml

**Lane 1** : SLAP2 transfected 293T cell lysate

**Lane 2** : Non-transfected 293T cell lysate

Lysates/proteins at 25 µg per lane.

### Secondary

**All lanes** : Goat Anti-Mouse IgG (H&L)-HRP at 1/2500 dilution

**Predicted band size:** 29 kDa

**Observed band size:** 29 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors