


Anti-ZKSCAN5 antibody ab66329

1 图像

概述

产品名称	Anti-ZKSCAN5抗体
描述	兔多克隆抗体to ZKSCAN5
经测试应用	适用于: WB, ELISA
种属反应性	与反应: Human 预测可用于: Horse, Guinea pig, Cat, Dog, Pig, Chimpanzee, African green monkey 
免疫原	A region within synthetic peptide: MIMTESREVI DLDPPAETSQ EQEDLFMKV EEEDCTWMQE YNPPTFETFY, corresponding to N terminal amino acids 1-50 of Human ZKSCAN5 Run BLAST with ExPASy Run BLAST with NCBI
阳性对照	Jurkat cell lysate.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None Constituents: 2% Sucrose, PBS
纯度	Protein A purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab66329** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab 评论	说明
WB		Use a concentration of 1 µg/ml. Detects a band of approximately 80 kDa (predicted molecular weight: 97 kDa). Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.


应用	Ab 评论	说明
----	-------	----

ELISA Use at an assay dependent concentration.
ELISA titre using peptide based assay: 1/62500.

靶标

功能	May be involved in transcriptional regulation.
序列相似性	Belongs to the krueppel C2H2-type zinc-finger protein family. Contains 13 C2H2-type zinc fingers. Contains 1 KRAB domain. Contains 1 SCAN box domain.
细胞定位	Nucleus.

Anti-ZKSCAN5 antibody 图像



Anti-ZKSCAN5 antibody (ab66329) at 1 µg/ml
(in 5% skim milk / PBS buffer) + Jurkat cell
lysate at 10 µg

Secondary

HRP conjugated anti-Rabbit IgG at 1/50000
dilution

Predicted band size : 97 kDa

Observed band size : 80 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors

