

## Product datasheet

# Anti-MTCP1 antibody ab67047

### 1 图像

#### 概述

产品名称	Anti-MTCP1抗体
描述	小鼠多克隆抗体to MTCP1
宿主	Mouse
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	Full length human MTCP1 protein
阳性对照	MTCP1 transfected 293T cell lysate.

#### 性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Constituent: PBS
纯度	Protein A purified
克隆	多克隆
同种型	IgG

#### 应用

Our [Abpromise guarantee](#) covers the use of **ab67047** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


应用	Ab评论	说明
WB		1/500 - 1/2500. Detects a band of approximately 12 kDa (predicted molecular weight: 7 kDa).

#### 靶标

功能	Enhances the phosphorylation and activation of AKT1 and AKT2.
----	---

<b>组织特异性</b>	Not found at a significant level in any tissue.
<b>疾病相关</b>	Note=Detected in T-cell leukemia bearing a t(X;14) translocation. Plays a key role in T-cell prolymphocytic leukemia.
<b>序列相似性</b>	Belongs to the TCL1 family.

## 图片



**All lanes** : Anti-MTCP1 antibody (ab67047)  
at 1/500 dilution

**Lane 1** : MTCP1 transfected 293T cell lysate

**Lane 2** : Non transfected 293T cell lysate

Lysates/proteins at 25 µg per lane.

### Secondary

**All lanes** : Goat Anti-Mouse IgG (H&L)-HRP  
Conjugate at 1/5000 dilution

**Predicted band size:** 7 kDa

**Observed band size:** 12 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors