

Product datasheet

Anti-SPAG4L antibody ab67057

1 图像

概述

产品名称	Anti-SPAG4L抗体
描述	小鼠多克隆抗体to SPAG4L
宿主	Mouse
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	SPAG4L (AAH26118, 1 a.a. ~ 380 a.a) full-length human protein.
阳性对照	Antibody reactive against mammalian transfected lysate.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None Constituents: Whole serum
纯度	Whole antiserum
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab67057** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/500 - 1/1000. Detects a band of approximately 43 kDa (predicted molecular weight: 43 kDa).

靶标

<b>相关性</b>	Sperm associated antigen 4 like protein. The function of this protein remains unknown.
<b>细胞定位</b>	Nuclear; integral component of membrane

## 图片



**All lanes** : Anti-SPAG4L antibody (ab67057)  
at 1/500 dilution

**Lane 1** : SPAG4L transfected 293T cell lysate

**Lane 2** : Non-transfected 293T cell lysate.

Lysates/proteins at 25 µg per lane.

### Secondary

**Lane 1** : Goat Anti-Mouse IgG HRP  
Conjugate at 1/5000 dilution

**Lane 2** : Goat Anti-Mouse IgG HRP  
Conjugate at 1/5000 dilution

**Predicted band size:** 43 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors