

Product datasheet

Anti-FOXS1 antibody ab67557

1 图像

概述

产品名称	Anti-FOXS1抗体
描述	小鼠多克隆抗体to FOXS1
宿主	Mouse
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	Full length FOXS1 protein (Human)
阳性对照	WB: FOXS1 transfected 293T cell lysate. ICC/IF: HeLa cell.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None Constituents: 1X PBS, pH 7.2
纯度	Protein A purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab67557** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		Use at an assay dependent dilution. Detects a band of approximately 36 kDa (predicted molecular weight: 36 kDa).

靶标

功能	Transcriptional repressor that suppresses transcription from the FASLG, FOXO3 and FOXO4 promoters. May have a role in the organization of the testicular vasculature.
序列相似性	Contains 1 fork-head DNA-binding domain.
细胞定位	Nucleus.

图片



All lanes : Anti-FOXS1 antibody (ab67557) at 1/500 dilution

Lane 1 : FOXS1 transfected 293T cell lysate

Lane 2 : Non transfected 293T cell lysate

Lysates/proteins at 25 µg per lane.

Secondary

All lanes : Goat Anti-Mouse IgG (H&L)-HRP Conjugate at 1/2500 dilution

Predicted band size: 36 kDa

Observed band size: 36 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors