

Product datasheet

Anti-CRTAC1 antibody ab67585

1 图像

概述

产品名称	Anti-CRTAC1抗体
描述	小鼠多克隆抗体to CRTAC1
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	Full length CRTAC1 protein (Human)
阳性对照	CRTAC1 transfected 293T cell lysate

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None Constituents: 1X PBS, pH 7.2
纯度	Protein A purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab67585** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/500 - 1/1000. Detects a band of approximately 50 kDa (predicted molecular weight: 71 kDa).

靶标

组织特异性	Expressed in the interterritorial matrix of articular deep zone cartilage (at protein level). Isoform 1
-------	---

and isoform 2 are expressed in brain. Isoform 1 is detected in lung and chondrocytes. Detected in cartilage, bone, cultured chondrocytes and lung, and at low levels in heart. Not detected in osteoblasts.

#### 序列相似性

Contains 1 EGF-like domain.

Contains 4 FG-GAP repeats.


#### 翻译后修饰

O-glycosylated.

#### 细胞定位

Secreted > extracellular space > extracellular matrix.

## 图片



**All lanes :** Anti-CRTAC1 antibody (ab67585)  
at 1/500 dilution

**Lane 1 :** CRTAC1 transfected 293T cell  
lysate

**Lane 2 :** Non-transfected 293T cell lysate

Lysates/proteins at 25 µg per lane.

#### Secondary

Goat Anti-Mouse IgG (H&L)-HRP Conjugate at  
1/2500 dilution

**Predicted band size :** 71 kDa

**Observed band size :** 50 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors