

Product datasheet

Anti-ZNF317 antibody ab67874

1 图像

概述

产品名称	Anti-ZNF317抗体
描述	小鼠多克隆抗体to ZNF317
宿主	Mouse
特异性	This antibody is reactive against mammalian transfected lysate.
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	ZNF317 (AAH78154.1, 1 a.a. ~ 595 a.a) full-length human protein.
阳性对照	ZNF317 transfected 293T cell lysate

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None Constituents: Whole serum
纯度	Whole antiserum
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab67874** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
----	------	----

WB

应用说明	WB: 1/500 - 1/1000. Detects a band of approximately 75 kDa (predicted molecular weight: 68 kDa).
------	--------------------------------------------------------------------------------------------------

Detection Antibody Pair.


Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

靶标

功能	May function as a transcription factor. May play an important role in erythroid maturation and lymphoid proliferation.
组织特异性	Isoform 1 and isoform 2 are ubiquitously expressed. Isoform 3 and isoform 4 are expressed only in lymphocytes, spleen and lung.
序列相似性	Belongs to the krueppel C2H2-type zinc-finger protein family. Contains 13 C2H2-type zinc fingers. Contains 1 KRAB domain.
细胞定位	Nucleus.

图片



Western blot - ZNF317 antibody (ab67874)

All lanes : Anti-ZNF317 antibody (ab67874)
at 1/500 dilution

Lane 1 : ZNF317 transfected 293T cell lysate

Lane 2 : Non-transfected 293T cell lysate

Lysates/proteins at 25 µg per lane.

Secondary

All lanes : Goat Anti-Mouse IgG (H&L)-HRP
Conjugate at 1/2500 dilution

Predicted band size: 68 kDa

Observed band size: 75 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you

- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors