

Product datasheet

Anti-PHF23 antibody ab67947

1 图像

概述

| | |
|-------|--|
| 产品名称 | Anti-PHF23抗体 |
| 描述 | 小鼠多克隆抗体to PHF23 |
| 宿主 | Mouse |
| 特异性 | Reactive against mammalian transfected lysate |
| 经测试应用 | 适用于: WB |
| 种属反应性 | 与反应: Human |
| 免疫原 | Full length protein, corresponding to amino acids 1-403 of Human PHF23 |
| 阳性对照 | Transfected 293T cell lysate |

性能

| | |
|------|---|
| 形式 | Liquid |
| 存放说明 | Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. |
| 存储溶液 | Preservative: None Constituents: Whole serum |
| 纯度 | Whole antiserum |
| 克隆 | 多克隆 |
| 同种型 | IgG |

应用

Our [Abpromise guarantee](#) covers the use of **ab67947** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| 应用 | Ab评论 | 说明 |
|----|------|--|
| WB | | 1/500 - 1/1000. Detects a band of approximately 70 kDa (predicted molecular weight: 44 kDa). |

图标

疾病相关

Note=A chromosomal aberration involving PHF23 is found in a patient with acute myeloid leukemia (AML). Translocation t(11;17)(p15;p13) with NUP98.


序列相似性

Belongs to the PHF23 family.
Contains 1 PHD-type zinc finger.

翻译后修饰

Phosphorylated upon DNA damage, probably by ATM or ATR.

图片



Western blot - PHF23 antibody (ab67947)

All lanes : Anti-PHF23 antibody (ab67947) at 1/500 dilution

Lane 1 : PHF23 transfected 293T lysate

Lane 2 : Non transfected 293T lysate

Lysates/proteins at 25 µg per lane.

Secondary

All lanes : Goat Anti-Mouse IgG (H&L)-HRP Conjugate at 1/2500 dilution

Predicted band size: 44 kDa

Observed band size: 70 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors