

Product datasheet

Anti-PHF23 antibody ab67947

1 图像

概述

产品名称	Anti-PHF23抗体
描述	小鼠多克隆抗体to PHF23
特异性	Reactive against mammalian transfected lysate
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	Full length protein, corresponding to amino acids 1-403 of Human PHF23
阳性对照	Transfected 293T cell lysate

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None Constituents: Whole serum
纯度	Whole antiserum
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab67947** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/500 - 1/1000. Detects a band of approximately 70 kDa (predicted molecular weight: 44 kDa).

靶标

## 疾病相关

Note=A chromosomal aberration involving PHF23 is found in a patient with acute myeloid leukemia (AML). Translocation t(11;17)(p15;p13) with NUP98.


## 序列相似性

Belongs to the PHF23 family.  
Contains 1 PHD-type zinc finger.

## 翻译后修饰

Phosphorylated upon DNA damage, probably by ATM or ATR.

## 图片



Western blot - PHF23 antibody (ab67947)

**All lanes** : Anti-PHF23 antibody (ab67947) at 1/500 dilution

**Lane 1** : PHF23 transfected 293T lysate

**Lane 2** : Non transfected 293T lysate

Lysates/proteins at 25 µg per lane.

### Secondary

Goat Anti-Mouse IgG (H&L)-HRP Conjugate at 1/2500 dilution

**Predicted band size** : 44 kDa

**Observed band size** : 70 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors