

Product datasheet

Anti-UEVLD antibody ab67948

1 图像

概述

产品名称	Anti-UEVLD抗体
描述	小鼠多克隆抗体to UEVLD
宿主	Mouse
特异性	Reactive against mammalian transfected lysate
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	Full length protein, corresponding to amino acids 1-357 of Human UEVLD
阳性对照	Transfected 293T cell lysate

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None Constituents: Whole serum
纯度	Whole antiserum
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab67948** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/500 - 1/1000. Detects a band of approximately 37 kDa (predicted molecular weight: 52 kDa).

图标

功能	Possible negative regulator of polyubiquitination.
组织特异性	Colon, colon carcinoma cell lines, normal cervical epithelium, carcinomas of the uterine cervix and peripheral blood leukocytes.
序列相似性	In the N-terminal section; belongs to the ubiquitin-conjugating enzyme family. UEV subfamily. In the C-terminal section; belongs to the LDH/MDH superfamily. Contains 1 UEV (ubiquitin E2 variant) domain.

图片



All lanes : Anti-UEVLD antibody (ab67948)
at 1/500 dilution

Lane 1 : UEVLD transfected 293T cell lysate

Lane 2 : Non transfected 293T cell lysate

Lysates/proteins at 25 µg per lane.

Secondary

All lanes : Goat Anti-Mouse IgG (H&L)-HRP
Conjugate at 1/2500 dilution

Predicted band size: 52 kDa

Observed band size: 37 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors