

Product datasheet

Anti-OTUB2 antibody ab68010

★★★★☆ 1 Abreviews 1 图像

概述

产品名称	Anti-OTUB2抗体
描述	小鼠多克隆抗体to OTUB2
特异性	This antibody is reactive against mammalian transfected lysate.
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	OTUB2 (NP_075601.1, 1 a.a. ~ 234 a.a) full-length human protein.
阳性对照	OTUB2 transfected 293T cell lysate

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None Constituents: Whole serum
纯度	Whole antiserum
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab68010** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB	★★★★☆	1/500 - 1/1000. Detects a band of approximately 25 kDa (predicted molecular weight: 27 kDa).

靶标

功能 Hydrolase that can remove conjugated ubiquitin from proteins in vitro and may therefore play an

important regulatory role at the level of protein turnover by preventing degradation. Mediates deubiquitination of both 'Lys-48'- and 'Lys-63'-linked polyubiquitin chains, with a preference for 'Lys-63'-linked polyubiquitin chains.

组织特异性


Widely expressed. Expressed at higher level in brain.

序列相似性

Belongs to the peptidase C65 family.

Contains 1 OTU domain.

Anti-OTUB2 antibody 图像



All lanes : Anti-OTUB2 antibody (ab68010) at 1/500 dilution

Lane 1 : OTUB2 transfected 293T cell lysate

Lane 2 : Non-transfected 293T cell lysate

Lysates/proteins at 25 µg per lane.

Secondary

Goat Anti-Mouse IgG (H&L)-HRP Conjugate at 1/2500 dilution

Predicted band size : 27 kDa

Observed band size : 25 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors