

Product datasheet

Anti-ZKSCAN1 antibody ab68254

1 图像

概述

产品名称	Anti-ZKSCAN1抗体
描述	小鼠多克隆抗体to ZKSCAN1
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	Full length ZKSCAN1 protein (Human)
阳性对照	ZKSCAN1 transfected 293T cell lysate

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None Constituents: PBS, pH 7.2
纯度	Protein A purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab68254** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/500 - 1/1000. Detects a band of approximately 50 kDa (predicted molecular weight: 64 kDa).

靶标

功能	May be involved in transcriptional regulation.
----	--

## 序列相似性

Belongs to the krueppel C2H2-type zinc-finger protein family.  
Contains 6 C2H2-type zinc fingers.  
Contains 1 KRAB domain.  
Contains 1 SCAN box domain.


## 翻译后修饰

Phosphorylated upon DNA damage, probably by ATM or ATR.

## 细胞定位

Nucleus.

## 图片



**All lanes :** Anti-ZKSCAN1 antibody  
(ab68254) at 1/500 dilution

**Lane 1 :** ZKSCAN1 transfected 293T cell  
lysate

**Lane 2 :** Non-transfected 293T cell lysate

Lysates/proteins at 25 µg per lane.

### Secondary

Goat Anti-Mouse IgG (H&L)-HRP Conjugate at  
1/2500 dilution

**Predicted band size :** 64 kDa

**Observed band size :** 50 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors