

Product datasheet

Anti-KLHDC4 antibody ab68350

1 图像

概述

产品名称	Anti-KLHDC4抗体
描述	小鼠单克隆抗体to KLHDC4
经测试应用	适用于: WB, ELISA
种属反应性	与反应: Human
免疫原	Recombinant full length KLHDC4 protein (Human), with a tag. MW of the tag alone is 26kDa.
阳性对照	KLHDC4 transfected 293T cell lysate

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None Constituents: 1X PBS, pH 7.2
纯度	Protein A purified
克隆	单克隆
同种型	IgG2a
轻链类型	kappa

应用


Our [Abpromise guarantee](#) covers the use of **ab68350** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/500 - 1/1000. Detects a band of approximately 70 kDa (predicted molecular weight: 58 kDa).
ELISA		Use at an assay dependent dilution.

靶标

## Anti-KLHDC4 antibody 图像



**All lanes :** Anti-KLHDC4 antibody (ab68350)  
at 1/500 dilution

**Lane 1 :** KLHDC4 transfected 293T cell  
lysate

**Lane 2 :** Non-transfected 293T cell lysate

Lysates/proteins at 25 µg/ml per lane.

### Secondary

Goat Anti-Mouse IgG (H&L)-HRP Conjugate at  
1/2500 dilution

**Predicted band size :** 58 kDa

**Observed band size :** 70 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors