

Product datasheet

Anti-STON1 antibody ab69590

1 图像

概述

产品名称	Anti-STON1抗体
描述	兔多克隆抗体to STON1
宿主	Rabbit
特异性	Detects endogenous levels of total STON1 protein.
经测试应用	适用于: WB, ELISA
种属反应性	与反应: Human
免疫原	Synthetic peptide from the C-terminal of human STON1.
阳性对照	K562 cell extract

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Preservative: 0.02% Sodium Azide Constituents: 50% Glycerol, PBS (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), 150mM Sodium chloride, pH 7.4
纯度	Immunogen affinity purified
纯化说明	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab69590** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/500 - 1/1000. Detects a band of approximately 83 kDa (predicted molecular weight: 83 kDa).


应用	Ab评论	说明
ELISA		1/10000.

## 靶标

**相关性** This protein may be involved in the endocytic machinery.

**细胞定位** Cell Membrane and Cytoplasmic

## 图片



**All lanes :** Anti-STON1 antibody (ab69590) at 1/500 dilution

**Lane 1 :** K562 cell extract

**Lane 2 :** K562 cell extract with peptide at 10 µg

Lysates/proteins at 30 µg per lane.

**Predicted band size:** 83 kDa

**Observed band size:** 83 kDa

**Additional bands at:** 34 kDa, 45 kDa, 55 kDa, 87 kDa. We are unsure as to the identity of these extra bands.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors