

Product datasheet

Anti-Keratin 36 antibody ab69832

1 图像

概述

产品名称	Anti-Keratin 36抗体
描述	小鼠多克隆抗体to Keratin 36
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	Keratin 36 (AAH43581, 1 a.a. ~ 417 a.a) human protein.
阳性对照	Keratin 36 transfected 293T cell lysate

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None Constituents: 1X PBS, pH 7.2
纯度	Protein G purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab69832** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
----	------	----

WB

应用说明	WB: Use at an assay dependent dilution. Detects a band of approximately 46 kDa (predicted molecular weight: 52 kDa). Not yet tested in other applications. Optimal dilutions/concentrations should be determined by the end user.
------	---

靶标


组织特异性

Expressed in the hair follicle.

序列相似性

Belongs to the intermediate filament family.

图片



All lanes : Anti-Keratin 36 antibody (ab69832) at 1/500 dilution

Lane 1 : Keratin 36 transfected 293T cell lysate

Lane 2 : Non transfected 293T cell lysate

Lysates/proteins at 25 µg per lane.

Secondary

Goat Anti-Mouse IgG (H&L)-HRP Conjugate at 1/2500 dilution

Predicted band size : 52 kDa

Observed band size : 46 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors