

Product datasheet

Anti-IFNA21 antibody ab70034

1 图像

概述

产品名称	Anti-IFNA21抗体
描述	小鼠多克隆抗体to IFNA21
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	Full length human IFNA21 protein (AAH69329.1).
阳性对照	IFNA21 transfected 293T cell lysate.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None Constituents: 1X PBS, pH 7.2
纯度	Protein G purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab70034** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/500 - 1/1000. Detects a band of approximately 22 kDa (predicted molecular weight: 22 kDa).

靶标


相关性 Produced by macrophages, IFN-alpha have antiviral activities. Interferon stimulates the

production of two enzymes: a protein kinase and an oligoadenylate synthetase.

细胞定位

Secreted

图片



All lanes : Anti-IFNA21 antibody (ab70034) at 1/500 dilution

Lane 1 : IFNA21 transfected 293T lysate

Lane 2 : Non-transfected 293T lysate

Lysates/proteins at 25 µg per lane.

Secondary

Goat Anti-Mouse IgG (H&L)-HRP at 1/2500 dilution

Predicted band size : 22 kDa

Observed band size : 22 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors