

Product datasheet

Goat Anti-Syrian Hamster IgG H&L (Alkaline Phosphatase) preadsorbed ab7147

概述

产品名称	山羊抗叙利亚仓鼠IgG H&L (Alkaline Phosphatase)预吸附二抗
宿主	Goat
靶标种属	Syrian hamster
经测试应用	适用于: ICC/IF, Dot blot, ELISA, IHC-P, IHC-Fr, Immunomicroscopy, WB
最小交叉反应	Chicken, Cow, Goat, Guinea pig, Horse, Human, Mouse, Rabbit, Rat, Sheep more details
免疫原	hamster IgG whole molecule
偶联物	Alkaline Phosphatase

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	pH: 7.60 Constituents: 0.88% Sodium chloride, 0.79% Tris HCl, 0.0095% Magnesium chloride, 0.0014% Zinc chloride, 50% Glycerol
纯度	IgG fraction
纯化说明	This product was prepared from monospecific antiserum by immunoaffinity chromatography using hamster IgG coupled to agarose beads.
共轭说明	Alkaline Phosphatase (Calf Intestine) (Molecular Weight 140,000 daltons)
克隆	多克隆

应用

Our [Abpromise guarantee](#) covers the use of **ab7147** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
ICC/IF		
Dot blot		
ELISA		
IHC-P		
IHC-Fr		
Immunomicroscopy		
WB		

应用说明 suitable for immunoblotting (western or dot blot), ELISA and immunohistochemistry as well as other phosphatase-antibody based enzymatic assays requiring lot-to-lot consistency

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors