

Product datasheet

Anti-RHBDL2 antibody ab71866

1 图像

概述

产品名称	Anti-RHBDL2抗体
描述	小鼠多克隆抗体to RHBDL2
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	Full length protein, corresponding to amino acids 1-121 of Human RHBDL2 (AAH13103.1, or Swissprot id Q9NX52-3 without the first 12 amino acids)
阳性对照	RHBDL2 transfected 293T lysate.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None Constituents: 1X PBS, pH 7.2
纯度	Protein A purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab71866** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


应用	Ab评论	说明
WB		1/500 - 1/1000. Detects a band of approximately 12 kDa (predicted molecular weight: 14 kDa).

靶标

**相关性** RHBDL2, rhomboid like protein 2, is involved in regulated intramembrane proteolysis and the subsequent release of functional polypeptides from their membrane anchors. There are 2 named isoforms produced by alternative splicing.

**细胞定位** Membrane; Multi pass membrane protein.

## 图片



**All lanes** : Anti-RHBDL2 antibody (ab71866)  
at 1/500 dilution

**Lane 1** : RHBDL2 transfected 293T lysate

**Lane 2** : Non transfected 293T lysate

Lysates/proteins at 25 µg per lane.

### Secondary

Goat Anti-Mouse IgG (H&L)-HRP at 1/2500  
dilution

**Predicted band size** : 14 kDa

**Observed band size** : 12 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors