

Product datasheet

Anti-RP11-484I6.3 antibody ab72928

1 图像

概述

产品名称	Anti-RP11-484I6.3抗体
描述	小鼠多克隆抗体to RP11-484I6.3
宿主	Mouse
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	Full length Human protein (NP_620134.2)
阳性对照	Transfected 293T cell lysate

性能

形式	Liquid
存放说明	Shipped at 4°C. After reconstitution store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Preservative: None Constituents: Whole serum
纯度	Whole antiserum
克隆	多克隆
同种型	IgG

应用


Our [Abpromise guarantee](#) covers the use of **ab72928** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/500 - 1/1000. Predicted molecular weight: 25 kDa.

靶标

## 图片



**All lanes :** Anti-RP11-484I6.3 antibody  
(ab72928) at 1/500 dilution

**Lane 1 :** RP11-484I6.3-transfected 293T cell  
lysate

**Lane 2 :** Non-transfected 293T cell lysate

Lysates/proteins at 25 µg per lane.

### Secondary

**All lanes :** HRP-conjugated goat anti-mouse  
at 1/2500 dilution

**Predicted band size:** 25 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors