

Product datasheet

Anti-GST antibody [4C10] ab73934

1 References

概述

| | |
|-------|--|
| 产品名称 | Anti-GST抗体[4C10] |
| 描述 | 小鼠单克隆抗体[4C10] to GST |
| 宿主 | Mouse |
| 经测试应用 | 适用于: WB, IP |
| 种属反应性 | N/A |
| 免疫原 | A GST-HS3 fusion protein. |
| 常规说明 | This antibody is effective for analyzing, characterizing and purifying fusion proteins containing the GST tag. |

性能

| | |
|--------------------|--|
| 形式 | Liquid |
| 存放说明 | Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. |
| 存储溶液 | Preservative: 0.03% Thimerosal (merthiolate) Constituents: PBS |
| 纯度 | Protein G purified |
| Primary antibody说明 | This antibody is effective for analyzing, characterizing and purifying fusion proteins containing the GST tag. |
| 克隆 | 单克隆 |
| 克隆编号 | 4C10 |
| 同种型 | IgG2b |

应用

Our [Abpromise guarantee](#) covers the use of **ab73934** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| 应用 | Ab评论 | 说明 |
|----|------|---|
| WB | | 1/1000. Predicted molecular weight: 26 kDa. |
| IP | | 1/150. |

靶标

相关性 GST (Glutathione S-Transferase) is a 26kDa protein encoded by the parasitic helminth *Schistosoma japonicum* and widely used in the pGEX family of GST plasmid expression vectors as a fusion protein with foreign proteins.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors