

Product datasheet

Anti-DGAT2L6 antibody ab74452

1 图像

概述

产品名称	Anti-DGAT2L6抗体
描述	兔多克隆抗体to DGAT2L6
宿主	Rabbit
特异性	ab74452 detects endogenous levels of total DGAT2L6 protein.
经测试应用	适用于: ELISA, WB
种属反应性	与反应: Human
免疫原	Synthetic peptide derived from an internal sequence of human DGAT2L6.
阳性对照	Jurkat cell extracts.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Preservative: 0.02% Sodium Azide Constituents: 50% Glycerol, PBS (without Mg ²⁺ and Ca ²⁺), 150mM Sodium chloride, pH 7.4
纯度	Immunogen affinity purified
纯化说明	ab74452 was affinity purified from rabbit antiserum by affinity chromatography using epitope specific immunogen.
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab74452** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


应用	Ab评论	说明
ELISA		1/10000.

应用	Ab评论	说明
WB		1/500 - 1/1000. Predicted molecular weight: 39 kDa.

靶标

功能	Probable acyltransferase uses fatty acyl-CoA as substrate (By similarity). Has no wax synthase activity to produce wax esters.
组织特异性	Expressed in all tissues tested except pancreas.
序列相似性	Belongs to the diacylglycerol acyltransferase family.
细胞定位	Endoplasmic reticulum membrane.

图片



All lanes : Anti-DGAT2L6 antibody (ab74452) at 1/500 dilution

Lane 1 : Jurkat cell extracts

Lane 2 : Jurkat cell extracts with immunising peptide at 10 µg

Lysates/proteins at 5 µg per lane.

Predicted band size: 39 kDa

Observed band size: 39 kDa

Additional bands at: 30 kDa, 80 kDa. We are unsure as to the identity of these extra bands.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise,

please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors