

Product datasheet

Anti-HOXA6 antibody [130C5a] ab74516

1 图像

概述

产品名称	Anti-HOXA6抗体[130C5a]
描述	小鼠单克隆抗体[130C5a] to HOXA6
宿主	Mouse
经测试应用	适用于: WB, Dot blot
种属反应性	与反应: Human
免疫原	Recombinant fragment from internal sequence of Human HOXA6.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Preservative: 0.05% Sodium Azide Constituents: 1% BSA, PBS, pH 7.4
纯度	Protein G purified
纯化说明	ab74516 was purified using protein G column chromatography from culture supernatant of hybridoma cultured in a medium containing bovine IgG-depleted (approximately 95%) fetal bovine serum and then filtered through a 0.22µm membrane.
克隆	单克隆
克隆编号	130C5a
同种型	IgG1

应用

Our [Abpromise guarantee](#) covers the use of **ab74516** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		Use at an assay dependent dilution. Predicted molecular weight: 26 kDa.
Dot blot		Use at an assay dependent dilution.

靶标

功能	Sequence-specific transcription factor which is part of a developmental regulatory system that provides cells with specific positional identities on the anterior-posterior axis.
序列相似性	Belongs to the Antp homeobox family. Contains 1 homeobox DNA-binding domain.
细胞定位	Nucleus.

图片



Western blot - Anti-HOXA6 antibody [130C5a]
(ab74516)

Staining of recombinant immunising protein using anti-HOXA6 antibody [130C5a] (ab74516).

Predicted band size: 26 kDa

Observed band size: 35 kDa

The molecular weight of the band on the western blot does not correspond to the molecular weight of the natural protein because only a fragment of the protein was used.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors