

Product datasheet

Anti-ENDOGL1 antibody ab74566

1 图像

概述

产品名称	Anti-ENDOGL1抗体
描述	兔多克隆抗体to ENDOGL1
经测试应用	适用于: WB, ELISA
种属反应性	与反应: Human
免疫原	Synthetic peptide derived from an internal sequence of human ENDOGL1.
阳性对照	Extracts from HUVEC, HT-29 and HepG2 cells.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Preservative: 0.02% Sodium Azide Constituents: 50% Glycerol, PBS (without Mg ²⁺ and Ca ²⁺), 150mM Sodium chloride, pH 7.4
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab74566** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/500 - 1/1000. Detects a band of approximately 41 kDa (predicted molecular weight: 41 kDa).
ELISA		1/10000.

靶标

功能 Endo/exonuclease with nicking activity towards supercoiled DNA, a preference for single

组织特异性

stranded DNA and 5'-3' exonuclease activity.

序列相似性


Ubiquitous.

细胞定位

Belongs to the DNA/RNA non-specific endonuclease family.

Mitochondrion inner membrane.

Anti-ENDOGL1 antibody 图像



All lanes : Anti-ENDOGL1 antibody (ab74566) at 1/500 dilution

Lane 1 : Extracts from HUVEC cells

Lane 2 : Extracts from HT-29 cells

Lane 3 : Extracts from HepG2 cells

Lane 4 : Extracts from HUVEC cells with immunising peptide at 5 µg

Lysates/proteins at 5 µg per lane.

Predicted band size : 41 kDa

Observed band size : 41 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors