


Product datasheet

# Anti-SNT2 antibody ab74861

1 图像

概述

产品名称	Anti-SNT2抗体
描述	兔多克隆抗体to SNT2
经测试应用	适用于: ELISA, WB
种属反应性	与反应: African green monkey 预测可用于: Human 
免疫原	Synthetic peptide derived from an internal sequence of human SNT2.
阳性对照	Extracts from COS-7 cells.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Preservative: 0.02% Sodium Azide Constituents: 50% Glycerol, PBS (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), 150mM Sodium chloride, pH 7.4
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab74861** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
ELISA		1/40000.
WB		1/500 - 1/1000. Predicted molecular weight: 54 kDa.

靶标

<b>相关性</b>	SNT2 is a substrate for the fibroblast growth factor receptor. It functions in linking FGF receptor stimulation to activators of Ras.
<b>细胞定位</b>	Plasma membrane

## 图片



Western blot - SNT2 antibody (ab74861)

**All lanes** : Anti-SNT2 antibody (ab74861) at 1/500 dilution

**Lane 1** : Extracts from COS-7 cells

**Lane 2** : Extracts from COS-7 cells with immunising peptide at 5 µg

Lysates/proteins at 5 µg per lane.

**Predicted band size** : 54 kDa

**Observed band size** : 54 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors