

Product datasheet

Anti-NOL10 antibody ab74986

1 图像

概述

| | |
|-------|---|
| 产品名称 | Anti-NOL10抗体 |
| 描述 | 兔多克隆抗体to NOL10 |
| 宿主 | Rabbit |
| 经测试应用 | 适用于: WB, ELISA |
| 种属反应性 | 与反应: Human |
| 免疫原 | A synthesized peptide derived from the C terminal of human NOL10. |
| 阳性对照 | Extracts from HUVEC cells. |

性能

| | |
|------|--|
| 形式 | Liquid |
| 存放说明 | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. |
| 存储溶液 | Preservative: 0.02% Sodium Azide Constituents: 50% Glycerol, PBS, 150mM Sodium chloride, pH 7.4 |
| 纯度 | Immunogen affinity purified |
| 克隆 | 多克隆 |
| 同种型 | IgG |

应用

Our [Abpromise guarantee](#) covers the use of **ab74986** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| 应用 | Ab评论 | 说明 |
|-------|------|--|
| WB | | 1/500 - 1/1000. Detects a band of approximately 80 kDa (predicted molecular weight: 80 kDa). |
| ELISA | | 1/40000. |

靶标

序列相似性


Belongs to the WD repeat NOL10/ENP2 family.

Contains 5 WD repeats.

细胞定位

Nucleus > nucleolus.

图片



All lanes : Anti-NOL10 antibody (ab74986) at 1/500 dilution

Lane 1 : Extracts from HUVEC cells

Lane 2 : Extracts from HUVEC cells plus 5 μ g immunizing peptide

Lysates/proteins at 5 μ g per lane.

Predicted band size: 80 kDa

Observed band size: 80 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors