


## Product datasheet

# Anti-KLKBL4 antibody ab76035

### 1 图像

#### 概述

产品名称	Anti-KLKBL4抗体
描述	兔多克隆抗体to KLKBL4
经测试应用	适用于: WB
种属反应性	与反应: Human 预测可用于: Horse 
免疫原	Synthetic peptide (Human)
阳性对照	Fetal muscle lysate

#### 性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None Constituents: 2% Sucrose, PBS
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

#### 应用

Our [Abpromise guarantee](#) covers the use of **ab76035** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 44 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

## 靶标

<b>序列相似性</b>	Belongs to the peptidase S1 family. Plasma kallikrein subfamily. Contains 1 peptidase S1 domain.
<b>结构域</b>	The serine protease domain lacks the essential His, Asp, and Ser residues of the catalytic triad at positions 83, 129 and 221 and is predicted to be inactive.
<b>细胞定位</b>	Secreted.

## 图片



Anti-KLKBL4 antibody (ab76035) at 1  $\mu$ g/ml  
(in 5% skim milk / PBS buffer) + fetal muscle  
lysate at 10  $\mu$ g

### Secondary

HRP-conjugated anti-Rabbit IgG at 1/50000  
dilution

**Predicted band size** : 44 kDa

**Observed band size** : 40 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors