

## Product datasheet

# Anti-LASS6 antibody ab79509

1 References 2 图像

### 概述

产品名称	Anti-LASS6抗体
描述	兔多克隆抗体to LASS6
宿主	Rabbit
特异性	Predicted not to cross-react with LASS5.
经测试应用	适用于: ICC/IF, WB, ELISA
种属反应性	与反应: Mouse, Rat, Human
免疫原	13 amino acid synthetic peptide from the N terminal region of the human LASS6
阳性对照	Rat brain tissue lysate

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Preservative: 0.02% Sodium Azide Constituents: PBS
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

### 应用

Our [Abpromise guarantee](#) covers the use of **ab79509** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
ICC/IF		Use a concentration of 20 µg/ml.
WB		Use a concentration of 1 - 2 µg/ml. Predicted molecular weight: 45 kDa.
ELISA		Use at an assay dependent dilution.

## 靶标

功能	May be involved in sphingolipid synthesis or its regulation.
序列相似性	Contains 1 homeobox DNA-binding domain. Contains 1 TLC (TRAM/LAG1/CLN8) domain.
翻译后修饰	Glycosylation on Asn-18 is not necessary for function.
细胞定位	Nucleus membrane. Endoplasmic reticulum membrane.

## 图片




**Lane 1** : Anti-LASS6 antibody (ab79509) at 1  $\mu\text{g/ml}$

**Lane 2** : Anti-LASS6 antibody (ab79509) at 2  $\mu\text{g/ml}$

**All lanes** : rat brain tissue lysate

Lysates/proteins at 15  $\mu\text{g}$  per lane.

**Predicted band size:** 45 kDa



Immunofluorescence of LASS6 in Mouse  
Brain cells using ab79509 at 5  $\mu\text{g/ml}$ .

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors