

Product datasheet

Anti-DRAK2 antibody ab806

1 图像

概述

产品名称	Anti-DRAK2抗体
描述	山羊多克隆抗体to DRAK2
宿主	Goat
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	Synthetic peptide: SLPNPHELVSDDL, corresponding to C terminal amino acids 359-371 of Human DRAK2. Run BLAST with ExPASy Run BLAST with NCBI
阳性对照	Jurkat cell lysates.
常规说明	DRAK2 is a serine threonine kinase with a death domain

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: 0.02% Sodium Azide Constituents: 0.5% BSA, 0.5mg/ml Tris, pH 7.3
纯度	IgG fraction
纯化说明	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Primary antibody说明	DRAK2 is a serine threonine kinase with a death domain
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab806** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		Use a concentration of 2 - 4 µg/ml. Can be blocked with Human DRAK2 peptide (ab22786) . Approx 45kDa band seen in Jurkat cell lysates.

靶标

功能	Acts as a positive regulator of apoptosis.
组织特异性	Highly expressed in placenta, lung, pancreas. Lower levels in heart, brain, liver, skeletal muscle and kidney.
序列相似性	Belongs to the protein kinase superfamily. CAMK Ser/Thr protein kinase family. DAP kinase subfamily. Contains 1 protein kinase domain.
翻译后修饰	Autophosphorylated.
细胞定位	Nucleus.

图片



staining (3µg/ml) of Jurkat lysate (RIPA buffer, 30µg total protein per lane). Primary incubated for 1 hour. Detected by western blot using chemiluminescence.

staining (3µg/ml) of Jurkat lysate (RIPA buffer, 30µg total protein per lane). Primary incubated for 1 hour. Detected by western blot using chemiluminescence.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise,

please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors