

Anti-Neisseria meningitidis antibody [5] ab80913

1 图像

概述

产品名称	Anti-Neisseria meningitidis抗体[5]
描述	小鼠单克隆抗体[5] to Neisseria meningitidis
宿主	Mouse
特异性	Recognizes recombinant Neisseria meningitidis LPDA protein
经测试应用	适用于: WB, ELISA
种属反应性	与反应: Recombinant fragment
免疫原	E.coli expressed recombinant protein corresponding to fragment (aa 332-469) Neisseria meningitidis LPDA
阳性对照	Recombinant Neisseria meningitidis LPDA protein

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: 0.02% Proclin 300 Constituents: Citrate Tris HCl, pH 7.0
纯度	Protein A purified
克隆	单克隆
克隆编号	5
同种型	IgG2b

应用

Our [Abpromise guarantee](#) covers the use of **ab80913** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		

应用	Ab评论	说明
----	------	----

ELISA

应用说明

ELISA: Use at an assay dependent dilution.
WB: 1/1000. Predicted molecular weight: 61 kDa.


This antibody has only been tested against the recombinant protein fragment.
Not yet tested in other applications.
Optimal dilutions/concentrations should be determined by the end user.

靶标

相关性

Neisseria meningitidis is a heterotrophic gram-negative diplococcal bacterium best known for its role in meningitis and other forms of meningococcal disease such as meningococemia. Strains of *Neisseria meningitidis* are characterised by the polysaccharide expressed on its capsule. Serogroups A, B and C cause 90% of meningococcal meningitis cases.

图片



Anti-*Neisseria meningitidis* antibody [5] (ab80913) at 1/1000 dilution + Recombinant *Neisseria meningitidis* LPDA

Predicted band size: 61 kDa
Observed band size: 40 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors