

Anti-Neisseria meningitidis antibody [5] ab80913

1 图像

概述

| | |
|-------|---|
| 产品名称 | Anti-Neisseria meningitidis抗体[5] |
| 描述 | 小鼠单克隆抗体[5] to Neisseria meningitidis |
| 特异性 | Recognizes recombinant Neisseria meningitidis LPDA protein |
| 经测试应用 | 适用于: WB, ELISA |
| 种属反应性 | 与反应: Recombinant fragment |
| 免疫原 | E.coli expressed recombinant protein corresponding to fragment (aa 332-469) Neisseria meningitidis LPDA |
| 阳性对照 | Recombinant Neisseria meningitidis LPDA protein |

性能

| | |
|------|---|
| 形式 | Liquid |
| 存放说明 | Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. |
| 存储溶液 | Preservative: 0.02% Proclin 300 Constituents: Citrate Tris HCl, pH 7.0 |
| 纯度 | Protein A purified |
| 克隆 | 单克隆 |
| 克隆编号 | 5 |
| 同种型 | IgG2b |

应用

Our [Abpromise guarantee](#) covers the use of **ab80913** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| 应用 | Ab评论 | 说明 |
|-------|------|----|
| WB | | |
| ELISA | | |

应用说明 ELISA: Use at an assay dependent dilution.

This antibody has only been tested against the recombinant protein fragment.

Not yet tested in other applications.


Optimal dilutions/concentrations should be determined by the end user.

靶标

相关性

Neisseria meningitidis is a heterotrophic gram-negative diplococcal bacterium best known for its role in meningitis and other forms of meningococcal disease such as meningococemia. Strains of *Neisseria meningitidis* are characterised by the polysaccharide expressed on its capsule. Serogroups A, B and C cause 90% of meningococcal meningitis cases.

Anti-*Neisseria meningitidis* antibody [5] 图像



Anti-*Neisseria meningitidis* antibody [5] (ab80913) at 1/1000 dilution + Recombinant *Neisseria meningitidis* LPDA

Predicted band size : 61 kDa

Observed band size : 40 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors