


Anti-SPOPL antibody ab81167

1 图像

概述

产品名称	Anti-SPOPL抗体
描述	兔多克隆抗体to SPOPL
经测试应用	适用于: ELISA, WB
种属反应性	与反应: Human 预测可用于: Mouse, Rat, Rabbit, Horse, Chicken, Guinea pig, Cow, Cat, Dog, Caenorhabditis elegans, Drosophila melanogaster, Zebrafish 
免疫原	Synthetic peptide: KGLDDESKDYLSLYLLLWSC PKSEVRKFK FSLLNAKREE TKAMESQRAY , corresponding to internal sequence amino acids 73-122 of Human SPOPL (NP_001001664). Run BLAST with ExPASy Run BLAST with NCBI
阳性对照	Jurkat cell lysate

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Preservative: None Constituents: 2% Sucrose, PBS
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab81167** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
ELISA		Use at an assay dependent concentration. ELISA titre using peptide based assay: 1/1562500.

应用	Ab评论	说明
----	------	----


WB Use a concentration of 1 µg/ml. Predicted molecular weight: 45 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

靶标

相关性 SPOPL belongs to the Tdpoz family. It contains 1 BTB (POZ) domain and 1 MATH domain. In complex with a cullin SPOPL may act in ubiquitination and proteasomal degradation processes.

细胞定位 Nuclear

Anti-SPOPL antibody 图像



Anti-SPOPL antibody (ab81167) at 1 µg/ml +
Jurkat cell lysate at 10 µg

Secondary

HRP conjugated anti-Rabbit IgG at 1/50000
dilution

Predicted band size : 45 kDa

Observed band size : 40 kDa

Additional bands at : 60 kDa. We are
unsure as to the identity of these extra bands.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors