


Product datasheet

# Anti-CTAGE5 antibody ab83248

1 图像

概述

<b>产品名称</b>	Anti-CTAGE5抗体
<b>描述</b>	兔多克隆抗体to CTAGE5
<b>经测试应用</b>	<b>适用于:</b> ELISA, WB
<b>种属反应性</b>	<b>与反应:</b> Human <b>预测可用于:</b> Mouse, Rat, Rabbit, Horse, Guinea pig, Cow, Cat, Dog, Pig, Chimpanzee 
<b>免疫原</b>	Synthetic peptide corresponding to a region within the internal sequence amino acids 539-588 (LLEGPLRLSP LLPGGGGRGS RGPGNPLDHQ ITNERGESSC DRLTDPHRAP) of Human CTAGE5 (NP_005921). <a href="#">Run BLAST with ExPASy</a> <a href="#">Run BLAST with NCBI</a>
<b>阳性对照</b>	293T cell lysate.

性能

<b>形式</b>	Liquid
<b>存放说明</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
<b>存储溶液</b>	Preservative: None Constituents: 2% Sucrose, PBS
<b>纯度</b>	Immunogen affinity purified
<b>克隆</b>	多克隆
<b>同种型</b>	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab83248** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


应用	Ab 评论	说明
ELISA		Use at an assay dependent concentration. ELISA titre using peptide based assay 1/62500.

应用	Ab 评论	说明
WB		Use a concentration of 1 µg/ml. Detects a band of approximately 91 kDa (predicted molecular weight: 91 kDa). Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

## 靶标

<b>功能</b>	Tumor-associated antigen.
<b>组织特异性</b>	Isoform 5A is expressed only in testis (at the protein level). Other isoforms are expressed in several other normal tissues, including brain, muscle and cranial skin. Isoform 5A (at protein level) and isoform 5B are expressed in cutaneous T-cell lymphoma (CTCL) cell lines, colorectal carcinomas, breast carcinomas and melanoma. Isoform 5B, but not isoform 5A, is expressed in head and neck squamous cell carcinoma.
<b>疾病相关</b>	Note=Autoantibodies against CTAGE5 are present in several cancer types, including benign meningioma and cutaneous T-cell lymphoma (CTCL).
<b>序列相似性</b>	Belongs to the cTAGE family.

## Anti-CTAGE5 antibody 图像



Anti-CTAGE5 antibody (ab83248) at 1 µg/ml  
(in 5% skim milk / PBS buffer) + 293T cell  
lysate at 10 µg

### Secondary

HRP conjugated anti-Rabbit IgG at 1/50000  
dilution

**Predicted band size :** 91 kDa

**Observed band size :** 91 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise,

please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors