

Product datasheet

Anti-C18orf25 antibody ab83567

4 图像

概述

<b>产品名称</b>	Anti-C18orf25抗体
<b>描述</b>	兔多克隆抗体to C18orf25
<b>宿主</b>	Rabbit
<b>经测试应用</b>	<b>适用于:</b> WB, ELISA
<b>种属反应性</b>	<b>与反应:</b> Human <b>预测可用于:</b> Mouse, Rat, Rabbit, Horse, Chicken, Cow, Cat, Dog, Pig 
<b>免疫原</b>	Synthetic peptide corresponding to a region within the internal sequence amino acids 216-265 (STSSDDDEE VSGSSKTITA EIPGHLDPGF LASDKTSAGN APLNEEINIA) of Human C18orf25 (NP_001008240) <a href="#">Run BLAST with ExPASy</a> <a href="#">Run BLAST with NCBI</a>
<b>阳性对照</b>	721_B cell lysate.

性能

<b>形式</b>	Liquid
<b>存放说明</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
<b>存储溶液</b>	Preservative: None Constituents: 2% Sucrose, PBS
<b>纯度</b>	Immunogen affinity purified
<b>克隆</b>	多克隆
<b>同种型</b>	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab83567** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明


应用	Ab评论	说明
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 43 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.
ELISA		Use at an assay dependent concentration. Titre using peptide based assay 1/62500.

## 靶标

### 翻译后修饰


Phosphorylated upon DNA damage, probably by ATM or ATR.

## 图片




Western blot analysis of Human fetal lung tissue lysate labeling C18orf25 with ab83567 at 1.0µg/ml.

Western blot - Anti-C18orf25 antibody (ab83567)




Western blot analysis of Human fetal liver tissue lysate labeling C18orf25 with ab83567 at 1.0µg/ml.

Western blot - Anti-C18orf25 antibody (ab83567)



Western blot analysis of Human fetal heart tissue lysate labeling C18orf25 with ab83567 at 1.0 µg/ml.



Anti-C18orf25 antibody (ab83567) at 1 µg/ml + 721\_B cell lysate at 10 µg

**Secondary**

HRP conjugated anti-Rabbit IgG at 1/50000 dilution

**Predicted band size:** 43 kDa

**Observed band size:** 39 kDa

Gel concentration 12%

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

**Our Abpromise to you: Quality guaranteed and expert technical support**

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

**Terms and conditions**

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors