

Product datasheet

Anti-DR3 antibody ab8413

1 图像

概述

产品名称	Anti-DR3抗体
描述	兔多克隆抗体to DR3
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Mouse, Rat, Human
免疫原	Synthetic peptide: QGGTRSPRCDGDFH , corresponding to amino acids 25/40 of Human DR3 (Peptide available as <a href="#">ab8445</a> .) <a href="#">Run BLAST with</a> <a href="#">Run BLAST with</a>
阳性对照	Jurkat total cell lysate A 59 kDa band should be detected.

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	PBS with 0.02% sodium azide
纯度	IgG fraction
克隆	多克隆
同种型	IgG
轻链类型	unknown

应用

Our [Abpromise guarantee](#) covers the use of **ab8413** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/500 - 1/1000. Detects a band of approximately 59 kDa.Can be blocked with <a href="#">Human DR3 peptide (ab8445)</a> .

## 靶标

<b>功能</b>	Receptor for TNFSF12/APO3L/TWEAK. Interacts directly with the adapter TRADD. Mediates activation of NF-kappa-B and induces apoptosis. May play a role in regulating lymphocyte homeostasis.
<b>组织特异性</b>	Abundantly expressed in thymocytes and lymphocytes. Detected in lymphocyte-rich tissues such as thymus, colon, intestine, and spleen. Also found in the prostate.
<b>序列相似性</b>	Contains 1 death domain. Contains 4 TNFR-Cys repeats.
<b>翻译后修饰</b>	Glycosylated.
<b>细胞定位</b>	Secreted and Cell membrane.

## 图片



Anti-DR3 antibody (ab8413) at 1/1000 dilution  
+ Jurkat total cell lysate

**Observed band size:** 59 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors