


Product datasheet

Anti-TAF5L antibody [935C5a] ab84871

1 图像

概述

产品名称	Anti-TAF5L抗体[935C5a]
描述	小鼠单克隆抗体[935C5a] to TAF5L
宿主	Mouse
经测试应用	适用于: WB, Dot blot
种属反应性	与反应: Recombinant fragment 预测可用于: Human 
免疫原	Recombinant fragment (tagged), corresponding to internal sequence amino acids 24-129 of Human TAF5L.
阳性对照	Immunising recombinant Human TAF5L protein.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: 0.05% Sodium Azide Constituents: 1% BSA, PBS, pH 7.4
纯度	Protein G purified
纯化说明	ab84871 was purified using protein G column chromatography from culture supernatant of hybridoma cultured in a medium containing bovine IgG-depleted (approximately 95%) fetal bovine serum. Filtered through a 0.22 µm membrane.
克隆	单克隆
克隆编号	935C5a
同种型	IgG2a

应用

Our [Abpromise guarantee](#) covers the use of **ab84871** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		Use at an assay dependent dilution. Predicted molecular weight: 66 kDa.
Dot blot		Use at an assay dependent dilution.

靶标

功能	Functions as a component of the PCAF complex. The PCAF complex is capable of efficiently acetylating histones in a nucleosomal context. The PCAF complex could be considered as the human version of the yeast SAGA complex.
序列相似性	Belongs to the WD repeat TAF5 family. Contains 6 WD repeats.
细胞定位	Nucleus.

图片



Anti-TAF5L antibody [935C5a] (ab84871) +
Immunogen (Recombinant Human TAF5L
fragment)

Predicted band size: 66 kDa

Observed band size: 36 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors