

## Product datasheet

# Anti-TAF5L antibody [935C5a] ab84871

### 1 图像

#### 概述

产品名称	Anti-TAF5L抗体[935C5a]
描述	小鼠单克隆抗体[935C5a] to TAF5L
宿主	Mouse
经测试应用	适用于: WB, Dot blot
种属反应性	与反应: Recombinant fragment 预测可用于: Human
免疫原	Recombinant fragment (tagged), corresponding to internal sequence amino acids 24-129 of Human TAF5L.
阳性对照	Immunising recombinant Human TAF5L protein.

#### 性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: 0.05% Sodium Azide Constituents: 1% BSA, PBS, pH 7.4
纯度	Protein G purified
纯化说明	ab84871 was purified using protein G column chromatography from culture supernatant of hybridoma cultured in a medium containing bovine IgG-depleted (approximately 95%) fetal bovine serum. Filtered through a 0.22 µm membrane.
克隆	单克隆
克隆编号	935C5a
同种型	IgG2a

#### 应用

Our [Abpromise guarantee](#) covers the use of **ab84871** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		Use at an assay dependent dilution. Predicted molecular weight: 66 kDa.
Dot blot		Use at an assay dependent dilution.

## 靶标

<b>功能</b>	Functions as a component of the PCAF complex. The PCAF complex is capable of efficiently acetylating histones in a nucleosomal context. The PCAF complex could be considered as the human version of the yeast SAGA complex.
<b>序列相似性</b>	Belongs to the WD repeat TAF5 family. Contains 6 WD repeats.
<b>细胞定位</b>	Nucleus.

## 图片



Anti-TAF5L antibody [935C5a] (ab84871) + Immunogen  
(Recombinant Human TAF5L fragment)

**Predicted band size:** 66 kDa

**Observed band size:** 36 kDa

[why is the actual band size different from the predicted?](#)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors