


Product datasheet

Anti-SCML4 antibody ab84929

1 图像

概述

| | |
|-------|---|
| 产品名称 | Anti-SCML4抗体 |
| 描述 | 兔多克隆抗体to SCML4 |
| 经测试应用 | 适用于: WB |
| 种属反应性 | 与反应: Human 预测可用于: Rat, Rabbit, Horse, Cow, Cat, Dog, Pig  |
| 免疫原 | Synthetic peptide, corresponding to a region within N terminal amino acids 71 - 120 (EPDLSSIPQDAATVPSLAAPQALTVCCLYINKQANAGPYLERKKVQQLPE H) of Human SCML4 (NP_932347) |
| 阳性对照 | OVCAR-3 cell lysate |

性能

| | |
|------|---|
| 形式 | Liquid |
| 存放说明 | Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. |
| 存储溶液 | Preservative: None Constituents: 2% Sucrose, PBS |
| 纯度 | Immunogen affinity purified |
| 克隆 | 多克隆 |
| 同种型 | IgG |

应用

Our [Abpromise guarantee](#) covers the use of **ab84929** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


| 应用 | Ab评论 | 说明 |
|----|------|--|
| WB | | Use a concentration of 1 µg/ml. Predicted molecular weight: 45 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T. |

靶标

相关性 SCML4 is a putative Polycomb group (PcG) protein. PcG proteins act by forming multiprotein complexes, which are required to maintain the transcriptionally repressive state of homeotic genes throughout development.

细胞定位 Nuclear

图片



Anti-SCML4 antibody (ab84929) at 1 µg/ml (in 5% skim milk / PBS buffer,) + OVCAR-3 cell lysate at 10 µg

Secondary

HRP conjugated anti-Rabbit IgG at 1/50000 dilution

Predicted band size : 45 kDa

Observed band size : 45 kDa

Additional bands at : 65 kDa. We are unsure as to the identity of these extra bands.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors