

Product datasheet

Anti-Melanoma antibody [HMB45 + M2-7C10 + M2-9E3 + T311]  
ab86735

1 References 1 图像

概述

产品名称	Anti-Melanoma抗体[HMB45 + M2-7C10 + M2-9E3 + T311]
描述	小鼠单克隆抗体[HMB45 + M2-7C10 + M2-9E3 + T311] to Melanoma
宿主	Mouse
经测试应用	适用于: IHC-P
种属反应性	与反应: Human
免疫原	Melanoma
阳性对照	Metastatic melanoma
常规说明	Isotype: IgG1/kappa + IgG2b + IgG2b + IgG2a

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Buffer with protein carrier and sodium azide (0.09%)
克隆	单克隆
克隆编号	HMB45 + M2-7C10 + M2-9E3 + T311

应用

Our [Abpromise guarantee](#) covers the use of **ab86735** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
----	------	----

IHC-P

应用说明	IHC-P: 1/50 - 1/100. Perform heat mediated antigen retrieval via the pressure cooker method before commencing with IHC staining protocol.
------	---

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

## 靶标


---

**相关性** Malignant melanoma is a malignant neoplasm of melanocytes, arising de novo or from a pre-existing benign nevus, which occurs most often in the skin but also may involve other sites.

**细胞定位** Membrane; Single-pass membrane protein

## 图片

---



Melanoma stained with Melanoma antibody cocktail ab86735 at 1/100 dilution.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Melanoma antibody [HMB45 + M2-7C10 + M2-9E3 + T311] (ab86735)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors